

ARG24070 Streptavidin (HRP) (Tag Free)

Package: 1 ml
Store at: -20°C

Summary


Product Description	HRP-conjugated Streptavidin, from Streptomyces avidinii.
Tested Application	ELISA, FLISA, WB
Target Name	Streptavidin
Expression System	Streptomyces avidinii
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	ELISA	1:4000 - 1:8000
	FLISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Purified
Buffer	PBS (pH 7.4) and 50% Glycerol.
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted protein at 4°C for up to 2-7 days. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. The protein solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG24070 Streptavidin (HRP) (Tag Free) WB image

Western blot: Rabbit IgG was resolved by electrophoresis under reducing conditions, transferred to PVDF membrane, and stained with [ARG21918](#) Goat anti-Rabbit IgG antibody (Biotin) (pre-adsorbed) followed by ARG24070 Streptavidin (HRP) (Tag Free) and chemiluminescent detection.



ARG24070 Streptavidin (HRP) (Tag Free) ELISA image

ELISA: The plate was coated with [ARG21905](#) Goat anti-Human IgG antibody (Biotin) (pre-adsorbed) and purified / unlabeled Rat IgG1 kappa. Biotin conjugated antibody and purified immunoglobulin were detected with serially diluted ARG24070 Streptavidin (HRP) (Tag Free).