

## ARG24075 Goat anti-Human IgG (Fc) antibody, pre-adsorbed

Package: 500 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes Human IgG (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, WB
Specificity	The antibody reacts with the Fc region of human IgG. The antibody is pre-adsorbed with Human IgG Fab, IgM and IgA, so the antibody may not react with Human IgG Fab, IgM and IgA, but may react with IgG from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Human IgG (Fc)
Species	Human
Immunogen	Human IgG
Conjugation	Un-conjugated

### Application Instructions

Pre Adsorbed	Human IgG Fab, IgM and IgA.	
Application table	Application	Dilution
	ELISA	1 µg/ml
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.