

ARG24089 anti-DDDDK tag antibody (PE)

Package: 100 µg
Store at: 4°C or -20°C

Summary

Product Description	PE-conjugated Rabbit polyclonal antibody recognizes DDDDK tag.
Tested Reactivity	Other
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	DDDDK tag
Species	Others
Immunogen	KLH conjugated synthetic peptide containing DYKDDDDK Tag.
Conjugation	PE

Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.4), 0.09% Sodium azide and 50% Glycerol.
Preservative	0.09% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody; Tag Internal Control antibody
---------------	--