

ARG24178 anti-Parainfluenza virus type 3 antibody

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody recognizes Parainfluenza virus type 3.
Tested Reactivity	Virus
Tested Application	ELISA, IHC-Fr
Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Target Name	Parainfluenza Virus type 3
Immunogen	Crude viral lysate, native.
Conjugation	Un-conjugated
Alternate Names	HPIV-3; HPIV3

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS and 0.09% Sodium Azide
Preservative	0.09% Sodium Azide
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Background	Human parainfluenza viruses (HPIVs) are the viruses that cause human parainfluenza. HPIVs are a paraphyletic group of four distinct single-stranded RNA viruses belonging to the Paramyxoviridae family. These viruses are closely associated with both human and veterinary disease. Virions are
------------	---

approximately 150–250 nm in size and contain negative sense RNA with a genome encompassing about 15,000 nucleotides.