

ARG30067 Free Prostate Specific Antigen ELISA Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10190	anti-PSA (free) antibody [108]	Mouse mAb	Hu	ELISA	100 µg
ARG10191	anti-PSA (free) antibody [103] (HRP)	Mouse mAb	Hu	ELISA	100 µl

Summary

Product Description	Free prostate specific antigen (fPSA) is the unbounded form of PSA. Study has shown that low percentage of fPSA in total PSA (tPSA) is associated with high risk of prostate cancer. Using percentage of fPSA in tPSA can enhance the predication of prostate cancer and avoid unnecessary biopsies in prostate screen. ARG30067 Free Prostate Specific Antigen ELISA Duos, includes a capture antibody, ARG10190 Free PSA antibody [108] and a HRP-conjugated tracer antibody, ARG10191 Free PSA antibody [103] (HRP), for studying Free prostate specific antigen (fPSA) protein expression level by ELISA.
Target Name	Free Prostate Specific Antigen
Alternate Names	Free Prostate Specific Antigen ELISA antibody; PSA (free) antibody; HRP-conjugated PSA (free) antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	ELISA Antibody Duo for Free Prostate Specific Antigen
Research Area	Cancer antibody