

ARG30279
IGFBP4 ELISA Antibody DuoPackage: 1 pair
Store at: -20°C**Component**

| Cat. No. | Component Name | Host clonality | Reactivity | Application | Package |
|----------|-------------------------------|----------------|------------|-------------|---------|
| ARG10611 | anti-IGFBP4 antibody [IBP3] | Mouse mAb | Hu | ELISA | 50 µg |
| ARG10612 | anti-IGFBP4 antibody [IBP163] | Mouse mAb | Hu | ELISA | 50 µg |

Summary

| | |
|-----------------|--|
| Target Name | IGFBP4 |
| Alternate Names | IGFBP4 ELISA antibody; IGFBP4 antibody |

Properties

| | |
|---------------------|--|
| Storage instruction | For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |