

**ARG30280**  
**IGFBP4 ELISA Antibody Duo**Package: 1 pair  
Store at: -20°C**Component**

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10613	anti-IGFBP4 antibody [IBP180]	Mouse mAb	Hu	ELISA	50 µg
ARG10614	anti-IGFBP4 antibody [IBP182]	Mouse mAb	Hu	ELISA	50 µg

**Summary**

Target Name	IGFBP4
Alternate Names	IGFBP4 ELISA antibody; IGFBP4 antibody

**Properties**

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.