

**ARG30281**  
**PLA2G7 / Lp-PLA2 ELISA Antibody Duo**Package: 1 pair  
Store at: -20°C

## Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10615	anti-PLA2G7 / Lp-PLA2 antibody [PL42cc]	Mouse mAb	Hu	ELISA, WB	50 µg
ARG10616	anti-PLA2G7 / Lp-PLA2 antibody [PL46cc]	Mouse mAb	Hu	ELISA, WB	50 µg

## Summary

Target Name	PLA2G7 / Lp-PLA2
Alternate Names	PLA2G7 / Lp-PLA2 ELISA antibody; PLA2G7 / Lp-PLA2 antibody

## Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.