

ARG30284
TIM1 / KIM1 ELISA Antibody DuoPackage: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10621	anti-TIM1 / KIM1 antibody [KIM70]	Mouse mAb	Hu	ELISA, WB	50 µg
ARG10622	anti-TIM1 / KIM1 antibody [KIM75]	Mouse mAb	Hu	ELISA, WB	50 µg

Summary

Target Name	TIM1 / KIM1
Alternate Names	TIM1 / KIM1 ELISA antibody; TIM1 / KIM1 antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.