

ARG30287 Influenza A Haemagglutinin H7 ELISA Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10628	anti-Influenza A Haemagglutinin H7 antibody [InA331]	Mouse mAb	Influenza A	ELISA, FLISA	50 µg
ARG10629	anti-Influenza A Haemagglutinin H7 antibody [InA334]	Mouse mAb	Influenza A	ELISA, FLISA, I-ELISA	50 µg


Summary

Target Name	Influenza A Haemagglutinin H7
Alternate Names	Influenza A Haemagglutinin H7 ELISA antibody; Influenza A Haemagglutinin H7 antibody

Properties


Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images




ARG10628 anti-Influenza A Haemagglutinin H7 antibody [InA331]
FLISA image

FLISA: Calibration curve for Influenza A haemagglutinin H7 in sandwich fluoroimmunoassay. Capture MAb: ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] - 1 µg/well. Detection MAb: (Eu-chelate labeled): ARG10628 anti-Influenza A Haemagglutinin H7 antibody [InA331] - 0.1 µg/well. Incubation time: 45 min. Antigens: Influenza A/Netherlands/219/03 H7N7 virus.



ARG10628 anti-Influenza A Haemagglutinin H7 antibody [InA331] FLISA image

FLISA: Calibration curve for Influenza A haemagglutinin H7 in sandwich fluoroimmunoassay. Capture MAb: ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] - 1 µg/well. Detection MAb: (Eu-chelate labeled): ARG10628 anti-Influenza A Haemagglutinin H7 antibody [InA331] - 0.1 µg/well. Incubation time: 45 min. Antigens: Recombinant H7 (A/Chicken/Netherlands/1/03) (aa 17-527).



ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] indirect ELISA image

Indirect ELISA: Titration curve of ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] specific to haemagglutinin H7. Antigens: Viruses - 200 ng/well and Recombinant H7 (A/Chicken/Netherlands/1/03)(aa 17-527) - 10 ng/well.