

ARG30290
Osteocalcin ELISA Antibody DuoPackage: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10528	anti-Osteocalcin antibody [2H9]	Mouse mAb	Hu	ELISA	125 µg
ARG10529	anti-Osteocalcin antibody [6F9]	Mouse mAb	Hu	ELISA	125 µg

Summary

Target Name	Osteocalcin
Alternate Names	Osteocalcin ELISA antibody; Osteocalcin antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.