

ARG40448 anti-MID1 antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes MID1
Tested Reactivity	Hu, Ms
Tested Application	ICC/IF, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	MID1
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 478-667 of Human MID1 (NP_000372.1).
Conjugation	Un-conjugated
Alternate Names	OSX; MIDIN; RING finger protein Midline-1; FXY; EC 6.3.2.-; E3 ubiquitin-protein ligase Midline-1; RNF59; Midin; TRIM18; Tripartite motif-containing protein 18; XPRF; ZNFXY; RING finger protein 59; GBBB1; OGS1; BBBG1; Putative transcription factor XPRF; OS

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Mouse heart	
Observed Size	85 kDa	


Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation


Gene Symbol	MID1
Gene Full Name	midline 1
Background	The protein encoded by this gene is a member of the tripartite motif (TRIM) family, also known as the 'RING-B box-coiled coil' (RBCC) subgroup of RING finger proteins. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein forms homodimers which associate with microtubules in the cytoplasm. The protein is likely involved in the formation of multiprotein structures acting as anchor points to microtubules. Mutations in this gene have been associated with the X-linked form of Opitz syndrome, which is characterized by midline abnormalities such as cleft lip, laryngeal cleft, heart defects, hypospadias, and agenesis of the corpus callosum. This gene was also the first example of a gene subject to X inactivation in human while escaping it in mouse. Multiple different transcript variants are generated by alternate splicing; however, the full-length nature of some of the variants has not been determined. [provided by RefSeq, Jul 2010]
Function	Has E3 ubiquitin ligase activity towards IGBP1, promoting its monoubiquitination, which results in deprotection of the catalytic subunit of protein phosphatase PP2A, and its subsequent degradation by polyubiquitination. [UniProt]
Calculated Mw	75 kDa
PTM	Phosphorylated on serine and threonine residues. [UniProt]
Cellular Localization	Cytoplasm. Cytoplasm, cytoskeleton. Cytoplasm, cytoskeleton, spindle. Note=Microtubule-associated. It is associated with microtubules throughout the cell cycle, co-localizing with cytoplasmic fibers in interphase and with the mitotic spindle and midbodies during mitosis and cytokinesis. [UniProt]

Images



ARG40448 anti-MID1 antibody ICC/IF image

Immunofluorescence: MCF7 cells stained with ARG40448 anti-MID1 antibody.



ARG40448 anti-MID1 antibody WB image

Western blot: 25 µg of Mouse heart lysate stained with ARG40448 anti-MID1 antibody at 1:1000 dilution.