

ARG42126
anti-SAP30BP antibodyPackage: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes SAP30BP
Tested Reactivity	Hu
Tested Application	ICC/IF, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	SAP30BP
Species	Human
Immunogen	KLH-conjugated synthetic peptide between aa. 160-188 of Human SAP30BP.
Conjugation	Un-conjugated
Alternate Names	HTRG; HCNGP; SAP30-binding protein; HTRP; Transcriptional regulator protein HCNGP

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:25
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Human placenta	
Observed Size	~ 46 kDa	


Properties

Form	Liquid
Purification	Purification with Protein A and immunogen peptide.
Buffer	PBS and 0.09% (W/V) Sodium azide.
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	SAP30BP
Gene Full Name	SAP30 binding protein
Function	Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes. [UniProt]
Calculated Mw	34 kDa
Cellular Localization	Nucleus. [UniProt]

Images



ARG42126 anti-SAP30BP antibody ICC/IF image

Immunofluorescence: HeLa cells stained with ARG42126 anti-SAP30BP antibody (green) at 1:25 dilution. Alexa Fluor® 555 Phalloidin (red) for cytoplasmic actin staining.



Human placenta

ARG42126 anti-SAP30BP antibody WB image

Western blot: 35 µg of Human placenta lysate stained with ARG42126 anti-SAP30BP antibody.