

ARG42525
anti-HA tag antibodyPackage: 100 µg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes HA tag
Tested Reactivity	Other
Tested Application	WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	HA tag
Species	Others
Immunogen	Purified recombinant peptide (MAYPYDVPDYANPMAYPYDVPDYANPMAYPYDVPDYA) derived from the Human Influenza Hemagglutinin protein.
Conjugation	Un-conjugated


Application Instructions

Application table	Application	Dilution
	WB	1:250 - 1:2000

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.05% Sodium azide and 20% Glycerol.
Preservative	0.05% Sodium azide
Stabilizer	20% Glycerol
Concentration	4 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



HA-BAG2
transfected 293

ARG42525 anti-HA tag antibody WB image

Western blot: HA-BAG2 transfected 293 cells. 40 μ g of cell lysate stained with ARG42525 anti-HA tag antibody at 1:1000 dilution.