

ARG42529 anti-Azurin antibody

Package: 100 µg
Store at: -20°C

Summary


Product Description	Goat Polyclonal antibody recognizes Azurin
Tested Reactivity	Bacteria
Tested Application	WB
Specificity	The antibody reacts with Burkholderia and Pseudomonas Azurin.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Azurin
Species	Bacteria
Immunogen	Recombinant Burkholderia Azurin protein.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:5000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.05% Sodium azide and 20% Glycerol.
Preservative	0.05% Sodium azide
Stabilizer	20% Glycerol
Concentration	3 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG42529 anti-Azurin antibody WB image

Western blot: 100 µg of total protein from Pseudomonas and 50 ng of His-Azurin recombinant protein. The blots were stained with ARG42529 anti-Azurin antibody at 1:2500 dilution.