

ARG42538 anti-Ty1 tag antibody

Package: 50 µg
Store at: -20°C

Summary


Product Description	Goat Polyclonal antibody recognizes Ty1 tag
Tested Reactivity	Other
Tested Application	WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Ty1 tag
Species	Others
Immunogen	Purified recombinant peptide (EVHTNQDPLD).
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:5000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.05% Sodium azide and 20% Glycerol.
Preservative	0.05% Sodium azide
Stabilizer	20% Glycerol
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG42538 anti-Ty1 tag antibody WB image

Western blot: 10 ng, 5 ng, 3 ng and 1 ng of NusA-Ty1 stained with ARG42538 anti-Ty1 tag antibody at 1:2500 dilution.