

ARG43147 anti-UBL3 antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes UBL3
Tested Reactivity	Hu, Ms, Rat
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	UBL3
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 1-117 of Human UBL3 (NP_009037.1).
Conjugation	Un-conjugated
Alternate Names	Membrane-anchored ubiquitin-fold protein; Ubiquitin-like protein 3; HsMUB; MUB; Protein HCG-1; HCG-1; PNSC1

Application Instructions

Application table	Application	Dilution
	IHC-P	1:50 - 1:200
	WB	1:1000 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Mouse brain	
Observed Size	~ 13 kDa	


Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	UBL3
Gene Full Name	ubiquitin-like 3
Calculated Mw	13 kDa
Cellular Localization	Cell membrane; Lipid-anchor. [UniProt]

Images



ARG43147 anti-UBL3 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Rat testis tissue stained with ARG43147 anti-UBL3 antibody at 1:100 dilution.



ARG43147 anti-UBL3 antibody WB image

Western blot: 25 µg of Mouse brain lysate stained with ARG43147 anti-UBL3 antibody at 1:1000 dilution.