

## ARG43260 anti-MEX3A antibody

Package: 50 µg  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes MEX3A
Tested Reactivity	Hu
Predict Reactivity	Ms, Rat
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	MEX3A
Species	Human
Immunogen	A 14-amino acid synthetic peptide within aa. 420-470 of Human MEX3A.
Conjugation	Un-conjugated
Alternate Names	RING finger and KH domain-containing protein 4; RKHD4; RNA-binding protein MEX3A; MEX-3A

### Application Instructions

Application table	Application	Dilution
	IHC-P	2 - 20 µg/ml
	WB	1 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	SK-N-SH	
Observed Size	~ 65 kDa	


### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.02% Sodium azide.
Preservative	0.02% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	MEX3A
Gene Full Name	mex-3 RNA binding family member A
Function	RNA binding protein, may be involved in post-transcriptional regulatory mechanisms. [UniProt]
Calculated Mw	54 kDa
PTM	Phosphorylated. [UniProt]
Cellular Localization	Cytoplasm. Nucleus. Cytoplasm, P-body. Note=Predominantly expressed in the cytoplasm and shuttles between the cytoplasm and the nucleus through the CRM1 export pathway. [UniProt]

Images




ARG43260 anti-MEX3A antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human brain tissue stained with ARG43260 anti-MEX3A antibody at 20 µg/ml dilution.



ARG43260 anti-MEX3A antibody WB image

Western blot: SK-N-SH cell lysate stained with ARG43260 anti-MEX3A antibody at 1 µg/ml dilution in the absence (left) and the presence of blocking peptide (right).



ARG43260 anti-MEX3A antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human brain tissue stained with ARG43260 anti-MEX3A antibody at 2.5 µg/ml dilution.