

ARG45000 anti-Lambda Light Chain antibody [RM110]

Package: 50 µg
Store at: -20°C

Summary

Product Description	ARG45000 Rabbit Monoclonal antibody [RM110] recognizes Lambda Light Chain.
Tested Reactivity	Ms
Tested Application	ELISA, WB
Specificity	This antibody reacts to the Lambda light chain of mouse immunoglobulins. No cross reactivity with the kappa light chain, human IgG, rat IgG, or goat IgG.
Host	Rabbit
Clonality	Monoclonal
Clone	RM110
Isotype	IgG
Target Name	Lambda Light Chain
Species	Mouse
Immunogen	Mouse IgM Lambda
Target Ig	IgM Lambda
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	0.005 - 0.2 µg/mL
	WB	0.1 - 0.5 µg/mL
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Buffer	PBS with 50% Glycerol, 1% BSA and 0.09% sodium azide
Preservative	0.09% sodium azide
Stabilizer	50% Glycerol, 1% BSA and 0.09%
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.