

## ARG45006 anti-IgG2b antibody [RM108]

Package: 50 µg  
Store at: -20°C

### Summary

Product Description	ARG45006 Rabbit Monoclonal [RM108] recognizes IgG2b.
Tested Reactivity	Ms
Tested Application	ELISA, WB
Specificity	This antibody reacts to the Fc region of mouse IgG2b. No cross reactivity with mouse IgG1, IgG2a, IgG3, IgM, IgA, IgE, human IgG, rat IgG, or goat IgG.
Host	Rabbit
Clonality	Monoclonal
Clone	RM108
Isotype	IgG
Target Name	IgG2b
Species	Mouse
Immunogen	Mouse IgG2b
Target Ig	IgG2b
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	0.005 - 0.2 µg/mL
	WB	0.1 - 0.5 µg/ml

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Buffer	PBS with 50% Glycerol, 1% BSA and 0.09% sodium azide
Preservative	0.09% sodium azide
Stabilizer	50% Glycerol, 1% BSA and 0.09%
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.