

ARG54493 anti-AU1 tag antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes AU1 tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, IHC-P, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	AU1 tag
Species	Others
Immunogen	AU1 (DTYRYI) coupled to KLH
Conjugation	Un-conjugated

Application Instructions

Application Note	ELISA: Coating: 2 - 10 µg/ml Primary: 0.03 - 1 µg/ml Western blot: Colorimetric: 0.1 - 1 µg/ml Chemiluminescent: 0.03 - 1 µg/ml IHC, ICC/IF: 0.5 - 5 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

Properties

Form	Liquid
Buffer	PBS (pH 7.2) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody
---------------	-------------------------------