

ARG54521 anti-Glu-Glu tag antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Glu-Glu tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Glu-Glu tag
Species	Others
Immunogen	glu-glu (EYMPME) coupled to KLH.
Conjugation	Un-conjugated

Application Instructions

Application Note	ELISA: Coating: 2 - 10 µg/ml Primary: 0.03 - 1 µg/ml Western blot: Colorimetric: 0.1 - 1 µg/ml Chemiluminescent: 0.03 - 1 µg/ml IHC, ICC/IF: 0.5 - 5 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

Properties

Form	Liquid
Buffer	PBS (pH 7.2) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody
---------------	-------------------------------