

## ARG54528 anti-GST tag antibody

Package: 100 µg  
Store at: -20°C

### Summary

|                     |                                                |
|---------------------|------------------------------------------------|
| Product Description | Rabbit Polyclonal antibody recognizes GST tag  |
| Tested Reactivity   | Other                                          |
| Tested Application  | ELISA, ICC/IF, IHC-P, WB                       |
| Host                | Rabbit                                         |
| Clonality           | Polyclonal                                     |
| Isotype             | IgG                                            |
| Target Name         | GST tag                                        |
| Species             | Others                                         |
| Immunogen           | GST from Schistosoma japonicum coupled to KLH. |
| Conjugation         | Un-conjugated                                  |

### Application Instructions

|                  |                                                                                                                                                                                                                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Note | ELISA: Coating: 2 - 10 µg/ml Primary: 0.03 - 1 µg/ml Western blot: Colorimetric: 0.1 - 1 µg/ml<br>Chemiluminescent: 0.03 - 1 µg/ml IHC, ICC/IF: 0.5 - 5 µg/ml<br>* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist. |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### Properties

|                     |                                                                                                                                                                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form                | Liquid                                                                                                                                                                                                                                                                                                                     |
| Buffer              | PBS (pH 7.2) and 0.1% Sodium azide                                                                                                                                                                                                                                                                                         |
| Preservative        | 0.1% Sodium azide                                                                                                                                                                                                                                                                                                          |
| Storage instruction | For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.                                                                                                                                                                                                                                                       |

### Bioinformation

|               |                                                    |
|---------------|----------------------------------------------------|
| Research Area | Controls and Markers antibody; Metabolism antibody |
|---------------|----------------------------------------------------|