

ARG54533 anti-HA tag antibody

Package: 100 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes HA tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, IP, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	HA tag
Species	Others
Immunogen	HA (YPYDVPDYA) conjugated to KLH
Conjugation	Un-conjugated

Application Instructions


Application Note	Western blot 0.03 - 1.0 µg/ml Immunoprecipitation 1 - 4 µg/mg lysate Immunocytochemistry 2.5 - 10 µg/ml ELISA Coating 2 - 10 µg/ml Primary 0.03 - 1 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

Properties

Form	Liquid
Purification	Antibody was purified by peptide affinity chromatography.
Buffer	PBS and 0.09% Sodium azide
Preservative	0.09% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody; Gene Regulation antibody
---------------	---



ARG54533 anti-HA tag antibody WB image

Western Blot: 200, 100, or 50ng of E. coli whole cell lysate expressing a multi-tag fusion protein stained with ARG54533 anti-HA tag antibody at 0.04ug/ml dilution.