

ARG54544
anti-T7 tag antibodyPackage: 50 µg
Store at: 4°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes T7 tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, IHC-P, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	T7 tag
Species	Others
Immunogen	T7 (MASMTGGQQMG) coupled to KLH.
Conjugation	HRP

Application Instructions

Application Note	ELISA: 0.01 - 0.1 µg/ml Western blot: Colorimetric: 0.1 - 1 µg/ml Chemiluminescent: 0.03 - 1 µg/ml IHC, ICC/IF: 0.5 - 2 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

Properties

Form	Liquid
Buffer	50mM HEPES (pH 7.1), 0.1M NaCl, 0.2% BSA and 0.01% Thimerosal
Preservative	0.01% Thimerosal
Stabilizer	0.2% BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody; Gene Regulation antibody
---------------	---