

ARG54882 anti-RUNX2 antibody

Package: 100 μl Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes RUNX2
Tested Reactivity	Hu
Tested Application	FACS, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	lgG
Target Name	RUNX2
Species	Human
Immunogen	KLH-conjugated synthetic peptide corresponding to aa. 445-474 surrounding S533 of Human RUNX2.
Conjugation	Un-conjugated
Alternate Names	OSF2; Core-binding factor subunit alpha-1; CLCD; AML3; Acute myeloid leukemia 3 protein; Polyomavirus enhancer-binding protein 2 alpha A subunit; CBFA1; SL3/AKV core-binding factor alpha A subunit; PEA2aA; Osteoblast-specific transcription factor 2; SL3-3 enhancer factor 1 alpha A subunit; PEBP2-alpha A; Oncogene AML-3; CBF-alpha-1; CCD; PEBP2aA; CCD1; OSF-2; PEA2-alpha A; Runt- related transcription factor 2

Application Instructions

Application table	Application	Dilution	
	FACS	1:10 - 1:50	
	IHC-P	Assay-dependent	
	WB	1:1000	
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.		
Positive Control	CEM		

Properties

Form	Liquid
Purification	Purification with Protein A and immunogen peptide.
Buffer	PBS and 0.09% (W/V) Sodium azide
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

R	io	inf	rm	ati	ion
				au	

Database links	GeneID: 860 Human
	Swiss-port # Q13950 Human
Gene Symbol	RUNX2
Gene Full Name	runt-related transcription factor 2
Background	This gene is a member of the RUNX family of transcription factors and encodes a nuclear protein with an Runt DNA-binding domain. This protein is essential for osteoblastic differentiation and skeletal morphogenesis and acts as a scaffold for nucleic acids and regulatory factors involved in skeletal gene expression. The protein can bind DNA both as a monomer or, with more affinity, as a subunit of a heterodimeric complex. Mutations in this gene have been associated with the bone development disorder cleidocranial dysplasia (CCD). Transcript variants that encode different protein isoforms result from the use of alternate promoters as well as alternate splicing. [provided by RefSeq, Jul 2008]
Function	Transcription factor involved in osteoblastic differentiation and skeletal morphogenesis. Essential for the maturation of osteoblasts and both intramembranous and endochondral ossification. CBF binds to the core site, 5'-PYGPYGGT-3', of a number of enhancers and promoters, including murine leukemia virus, polyomavirus enhancer, T-cell receptor enhancers, osteocalcin, osteopontin, bone sialoprotein, alpha 1(I) collagen, LCK, IL-3 and GM-CSF promoters. In osteoblasts, supports transcription activation: synergizes with SPEN/MINT to enhance FGFR2-mediated activation of the osteocalcin FGF-responsive element (OCFRE) (By similarity). Inhibits KAT6B-dependent transcriptional activation. [UniProt]
Research Area	Developmental Biology antibody; Gene Regulation antibody
Calculated Mw	57 kDa
PTM	Phosphorylated; probably by MAP kinases (MAPK). Phosphorylation by HIPK3 is required for the SPEN/MINT and FGF2 transactivation during osteoblastic differentiation (By similarity). Phosphorylation at Ser-451 by CDK1 promotes endothelial cell proliferation required for tumor angiogenesis probably by facilitating cell cycle progression. Isoform 3 is phosphorylated on Ser-340.
Cellular Localization	Nucleus.

Images



ARG54882 anti-RUNX2 antibody IHC-P image

Immunohistochemistry: Formalin-fixed and paraffin-embedded Human tonsil tissue stained with ARG54882 anti-RUNX2 antibody.



ARG54882 anti-RUNX2 antibody WB image

Western blot: 35 μg of CEM cell lysate stained with ARG54882 anti-RUNX2 antibody.



ARG54882 anti-RUNX2 antibody FACS image

Flow Cytometry: NCI-H460 cells stained with ARG54882 anti-RUNX2 antibody (right histogram) or without primary antibody control (left histogram), followed by incubation with FITC labelled secondary antibody.