

ARG55072 anti-CCDC90B antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes CCDC90B
Tested Reactivity	Hu
Predict Reactivity	Ms, Rat
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	CCDC90B
Species	Human
Immunogen	KLH-conjugated synthetic peptide corresponding to aa. 123-152 (Center) of Human CCDC90B.
Conjugation	Un-conjugated
Alternate Names	Coiled-coil domain-containing protein 90B, mitochondrial; MDS025; MDS011

Application Instructions

Application table	Application	Dilution
	IHC-P	Assay-dependent
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Human ovary	


Properties

Form	Liquid
Purification	Purification with Protein A and immunogen peptide.
Buffer	PBS and 0.09% (W/V) Sodium azide
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Database links	GeneID: 60492 Human Swiss-port # Q9GZT6 Human
Gene Symbol	CCDC90B
Gene Full Name	coiled-coil domain containing 90B
Research Area	Cell Biology and Cellular Response antibody; Metabolism antibody; Signaling Transduction antibody
Calculated Mw	30 kDa
Cellular Localization	Mitochondrion. Membrane; Single-pass membrane protein

Images



ARG55072 anti-CCDC90B antibody IHC-P image

Immunohistochemistry: Formalin-fixed and paraffin-embedded Human pancreas tissue stained with ARG55072 anti-CCDC90B antibody.



ARG55072 anti-CCDC90B antibody WB image

Western blot: 35 µg of Human ovary lysate stained with ARG55072 anti-CCDC90B antibody at 1:1000 dilution.