

ARG55161 anti-ADCK5 antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes ADCK5
Tested Reactivity	Hu
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	ADCK5
Species	Human
Immunogen	KLH-conjugated synthetic peptide corresponding to aa. 434-464 (C-terminus) of Human ADCK5.
Conjugation	Un-conjugated
Alternate Names	EC 2.7.11.-; Uncharacterized aarF domain-containing protein kinase 5

Application Instructions

Application table	Application	Dilution
	IHC-P	Assay-dependent
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	293	


Properties

Form	Liquid
Purification	Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Buffer	PBS and 0.09% (W/V) Sodium azide
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Database links	GeneID: 203054 Human Swiss-port # Q3MIX3 Human
Gene Symbol	ADCK5
Gene Full Name	aarF domain containing kinase 5
Function	The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr). [UniProt]
Research Area	Cell Biology and Cellular Response antibody
Calculated Mw	66 kDa
Cellular Localization	Membrane; Single-pass membrane protein

Images



ARG55161 anti-ADCK5 antibody IHC-P image

Immunohistochemistry: Formalin-fixed and paraffin-embedded Human cancer tissue stained with ARG55161 anti-ADCK5 antibody.



ARG55161 anti-ADCK5 antibody WB image

Western blot: 293 cell lysate stained with ARG55161 anti-ADCK5 antibody.