

## ARG55623 anti-RHOXF2B antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes RHOXF2B
Tested Reactivity	Hu
Tested Application	ICC/IF, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	RHOXF2B
Species	Human
Immunogen	KLH-conjugated synthetic peptide corresponding to aa. 95-121 (Center) of Human RHOXF2B.
Conjugation	Un-conjugated
Alternate Names	Rhox homeobox family member 2B

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:10 - 1:50
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	A549	


### Properties

Form	Liquid
Purification	Purification with Protein A and immunogen peptide.
Buffer	PBS and 0.09% (W/V) Sodium azide
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation


Database links	<a href="#">GeneID: 727940 Human</a> <a href="#">Swiss-port # POC7M4 Human</a>
Gene Symbol	RHOXF2B
Gene Full Name	Rhox homeobox family, member 2B
Research Area	Gene Regulation antibody
Calculated Mw	32 kDa
Cellular Localization	Nucleus {ECO:0000255 PROSITE- ProRule:PRU00108}

## Images



ARG55623 anti-RHOXF2B antibody ICC/IF image

Immunofluorescence: A549 cells were fixed with 4% PFA (20 min), permeabilized with Triton X-100 (0.1%, 10 min), then stained with ARG55623 anti-RHOXF2B antibody (green) at 1:25 dilution, 1 hour at 37°C. Cytoplasmic actin was counterstained with Alexa Fluor® 555 (red) conjugated Phalloidin (7 units/ml, 1 hour at 37°C).



ARG55623 anti-RHOXF2B antibody WB image

Western blot: 35 µg of A549 cell lysate stained with ARG55623 anti-RHOXF2B antibody.