

ARG55731 anti-V5 tag antibody

Package: 100 µg
Store at: -20°C

Summary


Product Description	Goat Polyclonal antibody recognizes V5 tag
Tested Reactivity	Other
Tested Application	WB
Specificity	Reacts specifically with V5-tagged recombinant fusion proteins expressed in transfected mammalian cells or produced by in vitro translation.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	V5 tag
Species	Others
Immunogen	Purified recombinant peptide around aa. 95-108 of the P/V proteins of the Paramyxovirus SV5
Conjugation	Un-conjugated
Alternate Names	V5 epitope; GKIPNPLLGLDST epitope

Application Instructions

Application table	Application	Dilution
	WB	1:500-1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.05% Sodium azide and 20% Glycerol
Preservative	0.05% Sodium azide
Stabilizer	20% Glycerol
Concentration	3 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG55731 anti-V5 tag antibody WB image

Western blot: 100 µg of MBP-V5 recombinant protein and 293 cells transfected with V5-mCherry stained with ARG55731 anti-V5 tag antibody at 1:1000 dilution.