

ARG56018 anti-CD57 antibody [HNK-1]

Package: 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [HNK-1] recognizes CD57
Tested Reactivity	Hu
Species Does Not React With	Rat
Tested Application	ICC/IF, IHC-P
Host	Mouse
Clonality	Monoclonal
Clone	HNK-1
Isotype	IgM, kappa
Target Name	CD57
Species	Human
Immunogen	Human peripheral blood mononuclear cells.
Conjugation	Un-conjugated
Alternate Names	Glucuronosyltransferase P; CD57; LEU7; GlcAT-P; GLCATP; HNK1; GlcUAT-P; NK1; Beta-1,3-glucuronyltransferase 1; EC 2.4.1.135; Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1; GLCUATP; UDP-GlcUA:glycoprotein beta-1,3-glucuronyltransferase; NK-1

Application Instructions

Application table	Application	Dilution
	ICC/IF	0.5 - 1 µg/ml
	IHC-P	2 - 4 µg/ml
Application Note	IHC-P: Antigen Retrieval: Boil tissue section in 10 mM Tris with 1 mM EDTA (pH 9.0) for 10-20 min, followed by cooling at RT for 20 min. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties


Form	Liquid
Purification	PEG precipitation
Buffer	PBS (pH 7.4), 0.05% Sodium azide and 0.1 mg/ml BSA
Preservative	0.05% Sodium azide
Stabilizer	0.1 mg/ml BSA
Concentration	0.2 mg/ml

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Database links	GeneID: 27087 Human Swiss-port # Q9P2W7 Human
Gene Symbol	B3GAT1
Gene Full Name	beta-1,3-glucuronyltransferase 1
Background	The protein encoded by this gene is a member of the glucuronyltransferase gene family. These enzymes exhibit strict acceptor specificity, recognizing nonreducing terminal sugars and their anomeric linkages. This gene product functions as the key enzyme in a glucuronyl transfer reaction during the biosynthesis of the carbohydrate epitope HNK-1 (human natural killer-1, also known as CD57 and LEU7). Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008]
Function	Involved in the biosynthesis of L2/HNK-1 carbohydrate epitope on glycoproteins. Can also play a role in glycosaminoglycan biosynthesis. Substrates include asialo-orosomucoid (ASOR), asialo-fetuin, and asialo-neural cell adhesion molecule. Requires sphingomyelin for activity: stearyl-sphingomyelin was the most effective, followed by palmitoyl-sphingomyelin and lignoceroyl-sphingomyelin. Activity was demonstrated only for sphingomyelin with a saturated fatty acid and not for that with an unsaturated fatty acid, regardless of the length of the acyl group (By similarity). [UniProt]
Calculated Mw	38 kDa
PTM	The soluble form derives from the membrane form by proteolytic processing.
Cellular Localization	Cell surface, cytoplasmic

Images



ARG56018 anti-CD57 antibody [HNK-1] IHC-P image

Immunohistochemistry: Formalin-fixed, paraffin-embedded Human tonsil stained with ARG56018 anti-CD57 antibody [HNK-1].