

ARG56112 anti-Melanoma Associated Antigen antibody [KBA.62]

Package: 50 µg

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [KBA.62] recognizes Melanoma Associated Antigen
Tested Reactivity	Hu
Tested Application	ICC/IF, IHC-P
Host	Mouse
Clonality	Monoclonal
Clone	KBA.62
Isotype	IgG1, kappa
Target Name	Melanoma Associated Antigen
Species	Human
Immunogen	Human KAL cells derived from lymph node metastasis of malignant melanoma.
Conjugation	Un-conjugated


Application Instructions

Application table	Application	Dilution
	ICC/IF	0.5 - 1 µg/ml
	IHC-P	0.5 - 1 µg/ml
Application Note	Antigen retrieval for IHC-P: Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties


Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS (pH 7.4), 0.05% Sodium azide and 0.1 mg/ml BSA
Preservative	0.05% Sodium azide
Stabilizer	0.1 mg/ml BSA
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Images



ARG56112 anti-Melanoma Associated Antigen antibody [KBA.62] IHC-P image

Immunohistochemistry: Formalin-fixed, paraffin-embedded Human melanoma stained with ARG56112 anti-Melanoma Associated Antigen antibody [KBA.62].



ARG56112 anti-Melanoma Associated Antigen antibody [KBA.62] IHC-P image

Immunohistochemistry: Formalin-fixed, paraffin-embedded Human tongue stained with ARG56112 anti-Melanoma Associated Antigen antibody [KBA.62].