

ARG56271
anti-Glycogen Synthase 1 antibodyPackage: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Glycogen Synthase 1
Tested Reactivity	Hu, Ms, Rat
Tested Application	ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Glycogen Synthase 1
Species	Human
Immunogen	Recombinant protein of Human Glycogen Synthase 1
Conjugation	Un-conjugated
Alternate Names	GSY; GYS; EC 2.4.1.11; Glycogen [starch] synthase, muscle

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	IHC-P	1:50 - 1:200
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	HeLa	


Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	GYS1
Gene Full Name	glycogen synthase 1 (muscle)
Background	The protein encoded by this gene catalyzes the addition of glucose monomers to the growing glycogen molecule through the formation of alpha-1,4-glycoside linkages. Mutations in this gene are associated with muscle glycogen storage disease. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Sep 2009]
Function	Transfers the glycosyl residue from UDP-Glc to the non-reducing end of alpha-1,4-glucan. [UniProt]
Calculated Mw	84 kDa
PTM	Phosphorylation at Ser-8 by AMPK inactivates the enzyme activity. Primed phosphorylation at Ser-657 (site 5) by CSNK2A1 and CSNK2A2 is required for inhibitory phosphorylation at Ser-641 (site 3a), Ser-645 (site 3b), Ser-649 (site 3c) and Ser-653 (site 4) by GSK3A and GSK3B (By similarity). Phosphorylated at Ser-641 by DYRK2, leading to inactivation (By similarity). Phosphorylated at Ser-641 by PASK, leading to inactivation; phosphorylation by PASK is inhibited by glycogen. Dephosphorylation at Ser-641 and Ser-645 by PP1 activates the enzyme.

Images



ARG56271 anti-Glycogen Synthase 1 antibody ICC/IF image

Immunofluorescence: HeLa cells stained with ARG56271 anti-Glycogen Synthase 1 antibody.




ARG56271 anti-Glycogen Synthase 1 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human stomach cancer stained with ARG56271 anti-Glycogen Synthase 1 antibody at 1:200 dilution.

ARG56271 anti-Glycogen Synthase 1 antibody WB image

Western blot: HeLa cell lysate stained with ARG56271 anti-Glycogen Synthase 1 antibody.



HeLa
