

ARG56452 anti-GST tag antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes GST tag
Tested Reactivity	Other
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	GST tag
Species	Others
Immunogen	Purified recombinant E. coli BL21(DE3) GST protein fragments expressed in E.coli.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	WB	1:10000


Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Concentration	0.7 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Calculated Mw	26 kDa
---------------	--------



ARG56452 anti-GST tag antibody WB image

Western blot: Purified recombinant Human GST-Lin28B expressed in E. coli, using ARG56452 anti-GST tag antibody at 1:10000 dilution.