

ARG56454 anti-GST tag antibody (HRP)

Package: 100 µl
Store at: -20°C

Summary

Product Description	HRP-conjugated Rabbit Polyclonal antibody recognizes GST tag
Tested Reactivity	Other
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	GST tag
Species	Others
Immunogen	Purified recombinant E. coli BL21(DE3) GST protein expressed in E.coli.
Conjugation	HRP

Application Instructions


Application table	Application	Dilution
	WB	1:5000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4), 0.1 mg/ml BSA and 50% Glycerol.
Stabilizer	0.1 mg/ml BSA and 50% Glycerol
Concentration	0.3 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Calculated Mw	26 kDa
---------------	--------



ARG56454 anti-GST tag antibody (HRP) WB image

Western blot: Recombinant GST fusion protein stained with ARG56454 anti-GST tag antibody (HRP) at 1:5000 dilution.