

ARG56545 anti-Nitrotyrosine antibody (Biotin)

Package: 50 µg
Store at: 4°C

Summary

Product Description	Mouse Monoclonal antibody recognizes Nitrotyrosine
Tested Reactivity	Other
Tested Application	ELISA, IHC-P, IP, WB
Specificity	This antibody has less than 5% reactivity with Chlorotyrosine.
Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Target Name	Nitrotyrosine
Species	Others
Immunogen	Peroxynitrite-treated KLH.
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	IHC-P	Assay-dependent
	IP	Assay-dependent
	WB	1 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Buffer	PBS (pH 7.2), 0.02% Sodium azide, 50% Glycerol and 0.1% BSA.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol and 0.1% BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.