

## ARG57943 anti-SKIV2L2 antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes SKIV2L2
Tested Reactivity	Hu, Ms, Rat
Tested Application	IHC-P, IP, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	SKIV2L2
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 1-304 of Human SKIV2L2 (NP_056175.3).
Conjugation	Un-conjugated
Alternate Names	ATP-dependent RNA helicase SKIV2L2; TRAMP-like complex helicase; fSAP118; KIAA0052; Mtr4; Superkiller viralicidic activity 2-like 2; Dob1; EC 3.6.4.13

### Application Instructions

Application table	Application	Dilution
	IHC-P	1:50 - 1:200
	IP	Assay-dependent
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Mouse spleen	
Observed Size	~ 125 kDa	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

## Bioinformation


Gene Symbol	SKIV2L2
Gene Full Name	superkiller viralicidic activity 2-like 2 (S. cerevisiae)
Function	May be involved in pre-mRNA splicing. Associated with the RNA exosome complex and involved in the 3'-processing of the 7S pre-rRNA to the mature 5.8S rRNA. [UniProt]
Calculated Mw	118 kDa
Cellular Localization	Nucleus, nucleolus, nucleoplasm. [UniProt]

## Images



ARG57943 anti-SKIV2L2 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human esophagus stained with ARG57943 anti-SKIV2L2 antibody at 1:100 dilution.



ARG57943 anti-SKIV2L2 antibody WB image

Western blot: 25 µg of Mouse spleen lysate stained with ARG57943 anti-SKIV2L2 antibody at 1:1000 dilution.