

## ARG57955 anti-ABI3BP antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes ABI3BP
Tested Reactivity	Hu
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	ABI3BP
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 29-280 of Human ABI3BP (NP_056244.2).
Conjugation	Un-conjugated
Alternate Names	NESHBP; Target of Nesh-SH3; ABI gene family member 3-binding protein; Nesh-binding protein; TARSH; NeshBP; Tarsh

### Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	HeLa	
Observed Size	~ 118 kDa	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	ABI3BP
Gene Full Name	ABI family, member 3 (NESH) binding protein
Calculated Mw	119 kDa
Cellular Localization	Secreted. [UniProt]

Images



HeLa

ARG57955 anti-ABI3BP antibody WB image

Western blot: 25 µg of HeLa cell lysate stained with ARG57955 anti-ABI3BP antibody at 1:1000 dilution.