

## ARG59129 anti-SFRS3 / SRSF3 antibody

Package: 50 µg  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes SFRS3 / SRSF3
Tested Reactivity	Hu, Ms, Rat
Predict Reactivity	Hm
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	SFRS3 / SRSF3
Species	Human
Immunogen	Synthetic peptide corresponding to aa. 1-29 of Human SFRS3 / SRSF3. (MHRDSCPLDCKVYVGNLGNNGNKTLEA).
Conjugation	Un-conjugated
Alternate Names	Splicing factor, arginine/serine-rich 3; Serine/arginine-rich splicing factor 3; SFRS3; Pre-mRNA-splicing factor SRP20; SRp20

### Application Instructions

Application table	Application	Dilution
	WB	0.1 - 0.5 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.05% Sodium azide and 5% BSA.
Preservative	0.05% Sodium azide
Stabilizer	5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	SRSF3
Gene Full Name	serine/arginine-rich splicing factor 3
Background	The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Two transcript variants, one protein-coding and the other non-coding, have been found for this gene. [provided by RefSeq, Sep 2010]
Function	May be involved in RNA processing in relation with cellular proliferation and/or maturation. May function as export adapter involved in mRNA nuclear export such as of histone H2A. Binds mRNA which is thought to be transferred to the NXF1-NXT1 heterodimer for export (TAP/NXF1 pathway); enhances NXF1-NXT1 RNA-binding activity. RNA-binding is semi-sequence specific. [UniProt]
Calculated Mw	19 kDa
PTM	Phosphorylated by CLK1, CLK2, CLK3 and CLK4. Extensively phosphorylated on serine residues in the RS domain. [UniProt]
Cellular Localization	Nucleus. Nucleus speckle. Cytoplasm. Note=Recruited to nuclear speckles following interaction with YTHDC1. [UniProt]

Images



ARG59129 anti-SFRS3 / SRSF3 antibody WB image

Western blot: 50 µg of Rat brain, 50 µg of Mouse brain and 40 µg of HeLa lysates stained with ARG59129 anti-SFRS3 / SRSF3 antibody at 0.5 µg/ml dilution.