

ARG59331 anti-STRAP / Unrip antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes STRAP / Unrip
Tested Reactivity	Hu, Ms, Rat
Predict Reactivity	Hm
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	STRAP / Unrip
Species	Human
Immunogen	Synthetic peptide corresponding to aa. 78-104 of Human STRAP / Unrip. (AADFTAKVWDAVSGDELMTLAHKHIVK)
Conjugation	Un-conjugated
Alternate Names	WD-40 repeat protein PT-WD; UNRIP; MAWD; UNR-interacting protein; Serine-threonine kinase receptor-associated protein; PT-WD; MAP activator with WD repeats

Application Instructions

Application table	Application	Dilution
	IHC-P	0.5 - 1 µg/ml
	WB	0.1 - 0.5 µg/ml
Application Note	IHC-P: Antigen Retrieval: By heat mediation. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.05% Sodium azide and 5% BSA.
Preservative	0.05% Sodium azide
Stabilizer	5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed

before use.


Note

For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	STRAP
Gene Full Name	serine/threonine kinase receptor associated protein
Function	The SMN complex plays a catalyst role in the assembly of small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Thereby, plays an important role in the splicing of cellular pre-mRNAs. Most spliceosomal snRNPs contain a common set of Sm proteins SNRPB, SNRPD1, SNRPD2, SNRPD3, SNRPE, SNRPF and SNRPG that assemble in a heptameric protein ring on the Sm site of the small nuclear RNA to form the core snRNP. In the cytosol, the Sm proteins SNRPD1, SNRPD2, SNRPE, SNRPF and SNRPG are trapped in an inactive 6S pICln-Sm complex by the chaperone CLNS1A that controls the assembly of the core snRNP. Dissociation by the SMN complex of CLNS1A from the trapped Sm proteins and their transfer to an SMN-Sm complex triggers the assembly of core snRNPs and their transport to the nucleus. STRAP plays a role in the cellular distribution of the SMN complex. Negatively regulates TGF-beta signaling but positively regulates the PDPK1 kinase activity by enhancing its autophosphorylation and by significantly reducing the association of PDPK1 with 14-3-3 protein. [UniProt]
Calculated Mw	38 kDa
Cellular Localization	Cytoplasm. Nucleus. Note=Localized predominantly in the cytoplasm but also found in the nucleus. [UniProt]

Images




ARG59331 anti-STRAP / Unrip antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human mammary cancer stained with ARG59331 anti-STRAP / Unrip antibody.



ARG59331 anti-STRAP / Unrip antibody WB image

Western blot: 50 µg of Rat brain, 40 µg of HepG2, 40 µg of SMMC-7721, 40 µg of U87, 50 µg of Human placenta and 40 µg of HeLa cell lysate stained with ARG59331 anti-STRAP / Unrip antibody at 0.5 µg/ml dilution.



ARG59331 anti-STRAP / Unrip antibody WB image

Western blot: 50 µg of sample under reducing conditions. HEK293, Jurkat, PC-3, THP-1, Rat brain, HEPA1-6 and Raw246.7 whole cell lysates stained with ARG59331 anti-STRAP / Unrip antibody at 0.5 µg/ml dilution, overnight at 4°C.