

ARG59457 anti-SNAI3 antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes SNAI3
Tested Reactivity	Hu
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	SNAI3
Species	Human
Immunogen	Synthetic peptide corresponding to aa. 155-178 of Human SNAI3. (FHCHKPYHTLAGLARHRQLHCHLQ)
Conjugation	Un-conjugated
Alternate Names	Zfp293; SMUC; Protein snail homolog 3; Zinc finger protein SNAI3; SNAIL3; Zinc finger protein 293; ZNF293

Application Instructions

Application table	Application	Dilution
	WB	0.1 - 0.5 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.05% Sodium azide and 5% BSA.
Preservative	0.05% Sodium azide
Stabilizer	5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	SNAI3
Gene Full Name	snail family zinc finger 3
Background	SNAI3 is a member of the SNAIL gene family, named for the Drosophila snail gene, which plays roles in mesodermal formation during embryogenesis (Katoh and Katoh, 2003 [PubMed 12579345]).[supplied by OMIM, Apr 2009]
Function	Seems to inhibit myoblast differentiation. Transcriptional repressor of E-box-dependent transactivation of downstream myogenic bHLHs genes. Binds preferentially to the canonical E-box sequences 5'-CAGGTG-3' and 5'-CACCTG-3' (By similarity). [UniProt]
Calculated Mw	32 kDa
Cellular Localization	Nucleus. [UniProt]

Images



ARG59457 anti-SNAI3 antibody WB image

Western blot: 40 µg of MCF-7 whole cell lysate stained with ARG59457 anti-SNAI3 antibody at 0.5 µg/ml dilution.