

**ARG59797**  
**anti-Telethonin antibody**Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes Telethonin
Tested Reactivity	Hu
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Telethonin
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 1-167 of Human Telethonin (NP_003664.1).
Conjugation	Un-conjugated
Alternate Names	telethonin; T-cap; CMH25; Telethonin; CMD1N; TELE; Titin cap protein; LGMD2G

### Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Rat heart and Mouse skeletal muscle	
Observed Size	19 kDa	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

Gene Symbol	TCAP
Gene Full Name	titin-cap
Background	Sarcomere assembly is regulated by the muscle protein titin. Titin is a giant elastic protein with kinase activity that extends half the length of a sarcomere. It serves as a scaffold to which myofibrils and other muscle related proteins are attached. This gene encodes a protein found in striated and cardiac muscle that binds to the titin Z1-Z2 domains and is a substrate of titin kinase, interactions thought to be critical to sarcomere assembly. Mutations in this gene are associated with limb-girdle muscular dystrophy type 2G. [provided by RefSeq, Jul 2008]
Function	Muscle assembly regulating factor. Mediates the antiparallel assembly of titin (TTN) molecules at the sarcomeric Z-disk. [UniProt]
Calculated Mw	19 kDa
Cellular Localization	Cytoplasm, myofibril, sarcomere. [UniProt]

## Images



ARG59797 anti-Telethonin antibody WB image

Western blot: 25 µg of Rat heart and Mouse skeletal muscle lysates stained with ARG59797 anti-Telethonin antibody at 1:1000 dilution.