

## ARG59868 anti-SFRS9 antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes SFRS9
Tested Reactivity	Hu, Ms
Tested Application	ICC/IF, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	SFRS9
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 1-221 of Human SFRS9 (NP_003760.1).
Conjugation	Un-conjugated
Alternate Names	Serine/arginine-rich splicing factor 9; Pre-mRNA-splicing factor SRp30c; Splicing factor, arginine/serine-rich 9; SRp30c; SFRS9

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Mouse testis and BT-474	
Observed Size	24 kDa	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	SRSF9
Gene Full Name	serine/arginine-rich splicing factor 9
Background	The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Two pseudogenes, one on chromosome 15 and the other on chromosome 21, have been found for this gene. [provided by RefSeq, Sep 2010]
Function	Plays a role in constitutive splicing and can modulate the selection of alternative splice sites. Represses the splicing of MAPT/Tau exon 10. [UniProt]
Calculated Mw	26 kDa
PTM	Extensively phosphorylated on serine residues in the RS domain. [UniProt]
Cellular Localization	Nucleus. Note=Cellular stresses such as heat shock may induce localization to discrete nuclear bodies termed SAM68 nuclear bodies (SNBs), HAP bodies, or stress bodies. Numerous splicing factors including SRSF1/SFRS1/SF2, SRSF7/SFRS7, SAFB and KHDRBS1/SAM68 accumulate at these structures, which may participate in the post-transcriptional regulation of mRNAs in stressed cells. [UniProt]

Images



ARG59868 anti-SFRS9 antibody ICC/IF image

Immunofluorescence: U2OS cells stained with ARG59868 anti-SFRS9 antibody at 1:100 dilution.



ARG59868 anti-SFRS9 antibody WB image

Western blot: 25 µg of Mouse testis and BT-474 cell lysates stained with ARG59868 anti-SFRS9 antibody at 1:1000 dilution.