


ARG62703
anti-CD106 / VCAM1 antibody [STA] (Biotin)

Package: 100 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Mouse Monoclonal antibody [STA] recognizes CD106 / VCAM1
Tested Reactivity	Hu
Tested Application	FACS
Specificity	

Images



ARG62703 anti-CD106 / VCAM1 antibody [STA] (Biotin) FACS image

Flow Cytometry: TNF alpha-stimulated HUVEC cells (red) and unstimulated HUVEC cells (black) stained with ARG62703 anti-CD106 / VCAM1 antibody [STA] (Biotin), followed by Streptavidin (APC).