

## ARG63228 anti-RAB11A / YL8 antibody

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes RAB11A / YL8
Tested Reactivity	Ms
Predict Reactivity	Hu, Rat, Cow
Tested Application	WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	RAB11A / YL8
Species	Human
Immunogen	TENKPKVQCCQNI
Conjugation	Un-conjugated
Alternate Names	Rab-11; Ras-related protein Rab-11A; YL8

### Application Instructions

Application table	Application	Dilution
	WB	1 - 3 µg/ml

**Application Note**  
WB: Recommend incubate at RT for 1h.  
\* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties


Form	Liquid
Purification	Purified from goat serum by antigen affinity chromatography.
Buffer	Tris saline (pH 7.3), 0.02% Sodium azide and 0.5% BSA.
Preservative	0.02% Sodium azide
Stabilizer	0.5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

---


Database links	<a href="#">GeneID: 53869 Mouse</a> <a href="#">Swiss-port # P62492 Mouse</a>
Background	The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2011]
Research Area	Cancer antibody; Signaling Transduction antibody
Calculated Mw	24 kDa

## Images



ARG63228 anti-RAB11A / YL8 antibody WB image

Western Blot: mouse heart lysate (RIPA buffer, 30µg total protein per lane) stained with ARG63228 anti-RAB11A / YL8 antibody at 2 µg/ml dilution.



ARG63228 anti-RAB11A / YL8 antibody WB image

Western Blot: 1). Mock transfection; 2) RAB11A(RC200352) expressing plasmid transfected HEK293 cell lysate stained with ARG63228 anti-RAB11A / YL8 antibody