

## ARG63281 anti-SNX26 / TCGAP antibody

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes SNX26 / TCGAP
Tested Reactivity	Hu
Predict Reactivity	Ms
Tested Application	IHC-P
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	SNX26 / TCGAP
Species	Human
Immunogen	SWSLHSEGGQTRSYC
Conjugation	Un-conjugated
Alternate Names	Tc10/CDC42 GTPase-activating protein; SNX26; Rho-type GTPase-activating protein 33; Sorting nexin-26; TCGAP; Rho GTPase-activating protein 33; NOMA-GAP

### Application Instructions

Application table	Application	Dilution
	IHC-P	2.5 µg/ml
Application Note	IHC-P: Antigen Retrieval: Steam tissue section in Citrate buffer (pH 6.0). * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified from goat serum by antigen affinity chromatography.
Buffer	Tris saline (pH 7.3), 0.02% Sodium azide and 0.5% BSA.
Preservative	0.02% Sodium azide
Stabilizer	0.5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

---

Database links [GeneID: 115703 Human](#)

[Swiss-port # O14559 Human](#)


**Background** This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. Alternative splice variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Feb 2010]

**Research Area** Signaling Transduction antibody

**Calculated Mw** 137 kDa

## Images

---



ARG63281 anti-SNX26 / TCGAP antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human brain tissue. Antigen Retrieval: Steam tissue section in Citrate buffer (pH 6.0). The tissue section was stained with ARG63281 anti-SNX26 / TCGAP antibody at 2.5 µg/ml dilution followed by AP-staining.