

ARG64401 anti-EBPL41L5 antibody

Package: 100 µg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes EBPL41L5
Tested Reactivity	Hu
Tested Application	FACS, ICC/IF, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	EBPL41L5
Species	Human
Immunogen	C-ENLPQSPGTDQHD
Conjugation	Un-conjugated
Alternate Names	Band 4.1-like protein 5; YMO1; YRT; BE37

Application Instructions

Application table	Application	Dilution
	FACS	10 µg/ml
	ICC/IF	10 µg/ml
	WB	1 - 3 µg/ml
Application Note	WB: Recommend incubate at RT for 1h. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Purified from goat serum by antigen affinity chromatography.
Buffer	Tris saline (pH 7.3), 0.02% Sodium azide and 0.5% BSA.
Preservative	0.02% Sodium azide
Stabilizer	0.5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation


Database links	GeneID: 57669 Human Swiss-port # Q9HCM4 Human
Gene Symbol	EPB41L5
Gene Full Name	erythrocyte membrane protein band 4.1 like 5
Function	May contribute to the correct positioning of tight junctions during the establishment of polarity in epithelial cells. [UniProt]
Research Area	Signaling Transduction antibody
Calculated Mw	82 kDa

Images




ARG64401 anti-EBPL41L5 antibody ICC/IF image

Immunofluorescence: A549 stained with ARG64401 anti-EBPL41L5 antibody at 10 µg/ml dilution.



ARG64401 anti-EBPL41L5 antibody WB image

Western Blot: Human Brain (Hippocampus) lysate (35 µg protein in RIPA buffer) stained with ARG64401 anti-EBPL41L5 antibody at 1 µg/ml dilution.



ARG64401 anti-EBPL41L5 antibody FACS image

Flow Cytometry: A549 stained with ARG64401 anti-EBPL41L5 antibody at 10ug/ml dilution.