

ARG65335 Mouse IgM Isotype Control antibody [PFR-03] (FITC)

Package: 100 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [PFR-03] as a negative control antibody for Mouse IgM
Species Does Not React With	Hu, Ms, Rat, Pig
Tested Application	FACS
Specificity	This mouse IgM monoclonal antibody (clone PFR-03) reacts with undefined epitope on a plant pathogen.
Host	Mouse
Clonality	Monoclonal
Clone	PFR-03
Isotype	IgM
Target Name	Mouse IgM
Conjugation	FITC

Application Instructions

Cross Reactivity Note Negative Species: Human, Pig, Mouse, Rat, Other species Not tested (expected negative)


Application table	Application	Dilution
	FACS	1 - 10 µg/ml

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification Note	The purified antibody is conjugated with Fluorescein isothiocyanate (FITC) under optimum conditions. The reagent is free of unconjugated FITC.
Buffer	TBS (pH 8.0), 15 mM Sodium azide and 0.2% BSA.
Preservative	15 mM Sodium azide
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG65335 Mouse IgM Isotype Control antibody [PFR-03] (FITC) FACS image

Flow Cytometry: Human peripheral blood stained with ARG65335 Mouse IgM Isotype Control antibody [PFR-03] (FITC) at 1 µg/ml dilution.