

ARG65349 Rat IgG Isotype Control antibody

Package: 100 µg
Store at: -20°C

Summary

Product Description	Rat Polyclonal antibody as a negative control antibody for Rat IgG
Tested Application	IHC-P, IP
Host	Rat
Clonality	Polyclonal
Isotype	IgG
Target Name	Rat IgG
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein A affinity column purified
Buffer	0.01M PBS (pH 7.4)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody
---------------	-------------------------------