

ARG65352 Donkey anti-Goat IgG antibody (HRP)

Package: 500 µl
Store at: -20°C

Summary

Product Description	HRP-conjugated Donkey Polyclonal antibody recognizes Goat IgG
Tested Reactivity	Goat
Tested Application	ELISA, IHC-P, WB
Specificity	Minor cross-reactivity to human in higher concentration and less than 1% cross-reactivity with mouse, rabbit and rat IgG is observed.
Host	Donkey
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG
Conjugation	HRP

Application Instructions

Application Note	IHC: 1:500-1:2500; ELISA: 1:5000-1:25000; WB: 1:5000-1:25000 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--


Properties

Form	Liquid
Buffer	PBS (pH 7.3), 0.01% Thimerosal, 50% Glycerol and 1% BSA.
Preservative	0.01% Thimerosal
Stabilizer	50% Glycerol and 1% BSA
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Keep the antibody in the dark and keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Goat IgG ELISA Kits; Related news: JunB as a novel gate keeper for prostate cancer progression 14-3-3n as a promising target for the treatment of Major Depression Disorder Has "Obesity gene" been found?
-----------	---


Images



ARG65352 Donkey anti-Goat IgG antibody (HRP) WB image

Western blot: AGS cells stained with [ARG56421 anti-COX2 antibody](#), [ARG63122 anti-ATF2 antibody](#) and ARG65352 Donkey anti-Goat IgG antibody (HRP).

From Anh Duy Do et al. Front Immunol. (2022), [doi: 10.3389/fimmu.2021.796177](#), Fig. 5.D.



ARG65352 Donkey anti-Goat IgG antibody (HRP) WB image

Western blot: Endometrial Stromal cells stained with [ARG63922 anti-EP4 prostaglandin Receptor antibody](#) at 1:1000 dilution and ARG65352 Donkey anti-Goat IgG antibody (HRP).

From Qingqing Huang et al. Reprod Med Biol. (2021), [doi: 10.1002/rmb2.12423](#), Fig. 3.C.